

MAIDESITE ELECTRIC NURSING BED

迈德斯特 MD-W01

Operation Mode: Electric

Size: 2080*960*520mm

Load Capacity $\leq 173\text{kg}$

Packing Size: 2130*1130*500mm

Backrest Lifting: $0\sim 60^\circ$

Leg Up & Down: $-55\sim 55^\circ$

Integral Forward & Backward Tilting: 15°

Integral Lifting Adjustment: 470~700mm

Upward and Downward

起背落腿

可满足饮水、就餐等需求

Provide drinking and eating support



AKL: 20902020975

Caring for Others

照顾别人