



EZ-GO PATIENT TRANSFER SLIDE EZ610



Protect Caregiver

Reduce lifting injuries by
62%



Comfort Patient

Reducing Pressure



Safe Money

Reducing cost
for Hospital

AKL: 10903020714

Product Name	EZ-GO Patient Transfer Slide
Model Number	EZ610
Place of Origin	Taiwan
Material	Polyester, PER
Size	140cm X 70cm
Color	Khaki
Product Weight	1 Kg
Features	<p>Multifunctional and foldable 6 way plus Patient Transfer Slide is an all around handy tool for caregivers</p> <ul style="list-style-type: none"> ■ Easy Way patient handling (9 handles). ■ 1 securing belt. ■ Anti slip mesh located at buttock area providing great stability during transfer. ■ Head support: Improve patient comfort and stability during transfer. ■ Multifunction: Transfer Slide + Assisted Beld + Soft Stretcher. ■ Can be used to assist emergency evacuation. ■ Heavy Duty: SGS tested, weight capacity up to 300 kg.