

LIST OF AIMTECH CURETTAGE INSTRUMENT SET

No.	Deskripsi	Tipe / Kode	Qty
1	Instrument tray, 310 x 210 x 10 mm	51-130-04	1
2	Sims uterine curette, sharp, rigid, Fig. 1, 7.5 mm, 26 cm	35-530-01	1
3	Sims uterine curette, sharp, rigid, Fig. 2, 8.5 mm, 26 cm	35-530-02	1
4	Sims uterine curette, sharp, rigid, Fig. 3, 10.5 mm, 26 cm	35-530-03	1
5	Sims uterine curette, sharp, rigid, Fig. 4, 12.0 mm, 26 cm	35-530-04	1
6	Sims uterine curette, sharp, rigid, Fig. 5, 14.0 mm, 26 cm	35-530-05	1
7	Sims uterine curette, sharp, rigid, Fig. 6, 15.0 mm, 26 cm	35-530-06	1
8	Sims scissors straight, blunt/blunt, 20 cm	11-120-20	1
9	Sims uterine sound, malleable, 32 cm	35-121-32	1
10	Gross-Maier sponge forceps, straight, 25 cm	15-702-25	1
11	Schroder uterine tenaculum forceps, 25 cm	35-204-25	2
12	Bozemann forceps, straight, 26 cm	15-716-26	2
13	Forster-Ballenger forceps, straight, 25 cm	15-736-25	1
14	Hegar uterine dilator set 8 pcs., sizes from 3.0/4.0 - 17.0/18.0 mm, in canvas roll	35-111-03	1
15	Mathieu foreign body forceps, WL 28 cm	36-080-28	1
16	Grave vaginal specula Fig. 1, 75 x 20 mm	35-013-01	5
17	Grave vaginal specula Fig. 2, 95 x 35 mm	35-013-02	5
18	Grave vaginal specula Fig. 3, 115 x 35 mm	35-013-03	1
19	Sims vaginal specula, Fig. 1, 65 x 25mm/70 x 30mm	35-019-01	1
20	Heywood-Smith forceps, straight, 20.5 cm	15-730-21	1
21	Standard dressing forceps, straight, 20 cm	13-000-20	1
22	Metal catheter, female, 10 Charr., 15 cm	36-110-10	1
23	Backhaus towel forceps, 13 cm	15-801-13	1
24	Surgical scissors, Fig. 2, straight, 14 cm, sharp/blunt	11-002-14	2
25	Surgical scissors, Fig. 2, curved, 14 cm, sharp/blunt	11-003-14	2
26	Littauer stitch scissors, straight, 14 cm	11-706-14	1
27	Lister bandage scissors, 14 cm	48-007-14	1
28	Braun-Stadler episiotomy scissors, 14 cm	11-139-14	1
29	Mayo-Hegar needle holder, delicate, 14 cm	17-004-14	1
30	Halsted-Mosquito forceps, straight, 1x2, 12 cm	15-160-12	2
31	Halsted-Mosquito forceps, curved, 1x2, 12 cm	15-161-12	2
32	Scalpel handle, standard No. 3	10-103-01	1
33	Surgical needle B, triangular, Fig. 10	17-803-10	2
34	Standard dressing forceps, straight, 14.5 cm	13-000-14	1
35	Dissecting forceps 1x2, 14.5 cm	13-100-14	1
36	Sims vaginal specula, Fig. 2, 70 x 30mm/70 x 35mm	35-019-02	1

37	Sims vaginal specula, Fig. 3, 85 x 35mm/90 x 40mm	35-019-03	1
----	---	-----------	---

Curettage Instrument Set



Instrument Tray

51-130-04 310 x 210 x 10 mm



Sims Uterine Curette

35-530-01 26,0cm, 10¼"



Sims Uterine Curette

35-530-02 26,0cm, 10¼"



Sims Uterine Curette

35-530-03 26,0cm, 10¼"



Sims Uterine Curette

35-530-04 26,0cm, 10¼"




Sims Uterine Curette

35-530-05 26,0cm, 10¼"



Sims Uterine Curette

35-530-06 26,0cm, 10¼"



Sims Scissors (Straight)

11-120-20 20,0 cm, 8"



Sims Uterine Sound

35-121-32 32,0cm, 12½"



Gross-Maier Forceps (Straight)

15-702-25 25,0 cm, 9¾"



Schröder Uterine Tenaculum Forceps


35-204-25 25,0 cm, 9¾"



Bozemann Forceps (Straight)

15-716-26 26,0cm, 10¼"

Curettage Instrument Set



**Förster-Ballenger
Forceps (Straight)**


15-736-25 25,0 cm, 9 3/4 "



Set 35-110-04 – 35-110-18 | 8 pcs.

Hegar Uterine Dilator

35-111-03 31/4 - 17/18 mm



Mathieu Foreign Body Forceps

36-080-28 28,0 cm, 11 "



lateral fixation screw

Grave Vaginal Specula

35-013-01 75 x 20 mm



lateral fixation screw

Grave Vaginal Specula

35-013-02 95 x 35 mm



lateral fixation screw


Grave Vaginal Specula

35-013-03 115 x 35 mm




Sims Vaginal Specula

35-019-01 70x30 | 65x25



**Heywood-Smith Forceps,
(Straight)**

15-730-21 20,5 cm, 8 "



**Dressing Forceps
(Straight)**


13-000-20 20,0 cm, 8 "



for women
10 Charr.


**Halsted-Mosquito Forceps
(Curved)**

36-110-10 15,0 cm, 6 "



Backhaus Towel Forceps


15-801-13 13,0 cm, 5 1/8 "



Surgical Scissors (Straight)

11-002-14 14,0 cm, 5 1/2 "

Curettage Instrument Set



Surgical Scissors (Curved)

11-003-14 14,0 cm, 5 ½ "




Littauer Stitch Scissors (Straight)

11-706-14 14,0 cm, 5 ½ "



Lister Bandage Scissors


48-007-14 14,0 cm, 5 ½ "



Perineum scissors
(Episiotomy)

Braun-Stadler Scissors

11-139-14 14,0 cm, 5 ½ "



delicate
pattern

Mayo-Hegar Needle Holder

17-004-14 14,0 cm, 5 ½ "



Halsted-Mosquito Forceps (Straight)

15-160-12 12,5 cm, 5 "




Halsted-Mosquito Forceps (Curved)

15-161-12 12,5 cm, 5 "



Scalpel Handle Std.3

10-103-01 12,0 cm, 4 ¾ "




Kreis curved 3/8 circle



Packing unit: 12
pieces


Surgical Needle "B" Triangular

17-803-10 Fig. 10



**Standard Dressing
Forceps (Straight)**

13-000-14 14,5 cm, 5 ¾ "



1 x 2

Dissecting Forceps

13-100-14 14,5 cm, 5 ¾ "

Curettage Instrument Set



Sims Vaginal Specula

35-019-02 70x35 | 70x30



Sims Vaginal Specula

35-019-03 90x40 | 85x35