

LIST OF AIMTECH ORTHOPEDIC INSTRUMENT SET

No.	Deskripsi	Tipe / Kode	Qty
1	Scalpel handle No. 3, solid, handle style as No. 4	10-103-15	1
2	Scalpel handle, standard No. 4	10-104-01	1
3	Mayo-Stille scissors, straight, 15 cm	11-170-15	1
4	Mayo-Stille scissors, curved, 15 cm	11-171-15	1
5	Surgical scissors, Fig. 1, curved, 14 cm, blunt/blunt	11-001-14	1
6	Metzenbaum scissors, curved, Fig. 2, 23 cm	11-203-23	1
7	Dissecting forceps 1x2, 14.5 cm	13-100-14	2
8	Dressing forceps delicate, straight, 14 cm	13-002-14	2
9	Dressing forceps delicate, straight, 18 cm	13-002-18	2
10	Adson dissecting, 1x2, forceps, 15 cm	13-142-15	2
11	Forster-Ballenger forceps, straight, 25 cm	15-736-25	1
12	Forster-Ballenger forceps, curved, 25 cm	15-737-25	2
13	Backhaus towel forceps, 11 cm	15-801-11	6
14	Kelly forceps, curved, 14 cm	15-053-14	2
15	Pean forceps, curved, 16 cm	15-001-16	2
16	Pean forceps, straight, 18 cm	15-000-18	2
17	Ochsner-Kocher forceps, straight, 16 cm	15-100-16	2
18	Ochsner-Kocher forceps, straight, 20 cm	15-100-20	2
19	Hegar needle holder, heavy pattern, 17.5 cm	17-000-18	2
20	Probe, d 2 mm, 14.5 cm	10-602-14	1
21	Universal wire scissors, curved, serrated 12 cm	11-721-12	1
22	Israel retractor, 4-prong, 40x40 mm 25.5 cm	16-017-04	2
23	Kelly Grip retractor, 65 x 50 mm, 26 cm	16-047-22	2
24	Richardson retractor, 52x22 mm, 24 cm	16-047-14	2
25	Langenbeck retractor, 55x20 mm, curved	16-020-02	1
26	Langenbeck retractor, 65x20 mm, curved	16-020-01	1
27	Volkman retractor, hollow handle, 6-prongs, blunt 22 cm	16-009-06	2
28	Volkman retractor, hollow handle, 4-prongs, blunt 22 cm	16-009-04	2
29	Volkman retractor, 4-prongs, sharp, 22 cm	16-001-04	2
30	Volkman retractor, 1-prong, sharp, 22 cm	16-001-01	2
31	Volkman retractor, hollow handle, 1-prong, blunt 22 cm	16-009-01	2
32	Weitlaner self Ret retractor, blunt, 3x4, 20 cm	16-244-20	2
33	Allis tissue forceps, 5x6, straight, 15 cm	15-766-15	2
34	Davis spatula, double, 24.5 cm	31-109-24	1
35	Langenbeck periosteal elevator, 16 mm, 19 cm	40-401-02	1
36	Volkman bone curette, Fig. 1, oval, 18.5 cm, Ferozell handle	40-479-05	1
37	Volkman bone curette, Fig. 2, oval, 18.5 cm, Ferozell handle	40-479-06	1

38	Volkman bone curette, Fig. 3, oval, 18.5 cm, Ferozell handle	40-479-07	1
39	Stille chisel, 10 mm, 20 cm	40-282-10	1
40	Stille chisel, 12 mm, 20 cm	40-282-12	1
41	Probe, d 2 mm, 16 cm	10-602-16	1
42	Doyen director, curved, 14 cm	10-613-14	1
43	Gross-Maier sponge forceps, curved, 20 cm	15-703-20	1
44	Poole suction cannula, d 10 mm, 22 cm	10-914-10	1
45	Redon guide needle curved, 12 Charr. with step formation	10-705-12	1
46	Redon guide needle curved, 14 Charr. with step formation	10-705-14	1
47	Round bowls, d 116 mm, H 50 mm, 0.35 l	51-103-11	2
48	Needle case perforated, with compartments 150 x 95 x 13 mm	56-112-02	1
49	Stille chisel, 15 mm, 20 cm	40-282-15	1
50	Stille chisel, 20 mm, 20 cm	40-282-20	1
51	Stille chisel, 25 mm, 20 cm	40-282-25	1
52	Gouge, 10 mm, 14 cm	40-264-10	1
53	Stille osteotome, curved, 12 mm, 20 cm	40-285-12	1
54	Stille osteotome, curved, 15 mm, 20 cm	40-285-15	1
55	Stille osteotome, curved, 20 mm, 20 cm	40-285-20	1
56	Stille osteotome, curved, 25 mm, 20 cm	40-285-25	1
57	Stille gouge, 10 mm, 20 cm	40-288-10	1
58	Stille gouge, 12 mm, 20 cm	40-288-12	1
59	Stille gouge, 15 mm, 20 cm	40-288-15	1
60	Stille gouge, 20 mm, 20 cm	40-288-20	1
61	Stille gouge, 25 mm, 20 cm	40-288-25	1
62	Mallet with plastic handle, d 42 / 30 mm 440 g, 26 cm	40-228-42	1
63	Hackstock flat nose pliers, with groove, 16 cm	40-852-16	1
64	Liston bone cutting forceps, straight, 14 cm	40-706-14	1
65	Liston bone cutting forceps, angled, 20 cm	40-707-19	1

Orthopedic Instrument Set




Scalpel Handle Std 3 (Solid)

10-103-15 12,5 cm, 8 "




Scalpel Handle Std.3

10-104-01 13,5 cm, 5 ¼ "




Mayo-Stillie Scissors

11-170-15 15,0 cm, 6 "




Surgical Scissors (Curved)

11-001-14 14,0 cm, 5 ½ "




Metzenbaum Scissors (Curved)

11-203-23 23,0 cm, 8 "




Dissecting Forceps

13-100-14 14,5 cm, 5 ¾ "



Dressing Forceps

13-002-14 14,5 cm, 5 ¾ "




Dressing Forceps

13-002-18 18,0 cm, 7 "




Adson Dissecting Forceps

13-142-15 15,0 cm, 6 "




Förster-Ballenger Forceps (Straight)

15-736-25 25,0 cm, 9 ¾ "



Förster-Ballenger Forceps (Curved)


15-737-25 25,0 cm, 9 ¾ "



Backhaus Towel Forceps


15-801-11 11,0 cm, 4 ⅜ "

Orthopedic Instrument Set




Kelly Forceps (Curved)

15-053-14 14,0 cm, 5 1/2 "



Pean Forceps (Curved)

15-001-16 16,0 cm, 6 1/4 "




Pean Forceps (Straight)

15-000-18 18,0 cm, 7 "



Ochsner-Kocher Forceps (Straight)

15-100-16 16,0 cm, 6 1/4 "



Ochsner-Kocher Forceps (Straight)

15-100-20 20,0 cm, 8 "




Hegar Needle Holder

17-000-18 17,5 cm, 7 "




Probe, Ø 2 mm

10-602-14 14,5 cm, 5 3/4 "



Universal Wire Scissors (Curved)

11-721-12 12,0 cm, 4 3/4 "



Israel Retractor, 4 Prongs

16-017-04 25,5 cm, 10 "




Kelly Grip Retractor

16-047-24 26,0cm, 10 1/4 "



Richardson Retractor

16-047-14 24,0 cm, 9 1/2 "




Langenbeck Retractor

16-020-02 22,0 cm, 8 3/4 "


Orthopedic Instrument Set

65 x 20 mm




Langenbeck Retractor

16-020-01 21,0 cm, 8 1/4"




**Volkman Retractor
6 Prongs**

16-009-06 22,0 cm, 8 3/4"




**Volkman Retractor
4 Prongs**

16-009-04 22,0 cm, 8 3/4"




**Volkman Retractor
4 Prongs**

16-001-04 22,0 cm, 8 3/4"




**Volkman Retractor
1 Prong**

16-001-01 22,0 cm, 8 3/4"



**Volkman Retractor
1 Prong**

16-009-01 22,0 cm, 8 3/4"



**Weitlaner Self Ret
Retractor,**

16-244-20 20,0 cm, 8"



Allis Tissue Forceps

15-766-15 15,0 cm, 6"



Davis Spatula (Double)

31-109-24 24,5 cm, 9 2/3"



Langenbeck Periosteal Elevator Volkman Bone Currette

40-401-02 19,0 cm, 7 1/2"



40-479-05 18,5 cm, 7 1/4"



Volkman Bone Currette

40-479-06 18,5 cm, 7 1/4"

Orthopedic Instrument Set




Fig. 3



Volkman Bone Curette

40-479-07 18,5 cm, 7 1/4"



10,0 mm

Stille Chisel

40-282-10 20,0 cm, 8"



12,0 mm

Stille Chisel

40-282-12 20,0 cm, 8"



Ø 2,0 mm

Probe, Ø 2 mm

10-602-16 16,0 cm, 6 1/4"



Doyen Director (Curved)

10-613-14 14,5 cm, 5 3/4"



Gross-Maier Forceps (Curved)

15-703-20 20,0 cm, 8"



1,3 mm




10 mm

with tube to be
screwed off

Pooler Suction Canula


10-914-10 22,0 cm, 8 5/8"



12 Charr.
Knife edge and step
formation

**Redon Guide Needle
(Curved)**

10-705-12 16,0 cm, 6 1/4"



14 Charr.
Knife edge and step
formation

**Redon Guide Needle
(Curved)**

10-705-14 16,0 cm, 6 1/4"



Round Bowls

51-103-11; Ø 116 x 50 mm; 0,35 l



**Needle Case
Compartments**

56-112-02 150 x 95 x 13 mm




15,0 mm

Stille Chisel


40-282-15 20,0 cm, 8"

Orthopedic Instrument Set




Stille Chisel

40-282-20 20,0 cm, 8 "




Stille Chisel

40-282-25 20,0 cm, 8 "




Gouge

40-264-10 14,0 cm, 5 ½ "



Stille Osteotome (Curved)

40-285-12 20,0 cm, 8 "



Stille Osteotome (Curved)

40-285-15 20,0 cm, 8 "




Stille Osteotome (Curved)

40-285-20 20,0 cm, 8 "




Stille Osteotome (Curved)

40-285-25 20,0 cm, 8 "




Stille Gouge

40-288-10 20,0 cm, 8 "




Stille Gouge

40-288-12 20,0 cm, 8 "




Stille Gouge

40-288-15 20,0 cm, 8 "



Stille Gouge

40-288-20 20,0 cm, 8 "



Stille Gouge

40-288-25 20,0 cm, 8 "

Orthopedic Instrument Set



HW 440 g
TW 640 g
Ø 30 | 42 mm



Mallet with Plastic Handle

40-228-42 26,0cm, 10 1/4"



Hackstock Flat Nose Pliers with Groove

40-852-16 16,0 cm, 6 1/4"



Liston Bone Cutting Forceps (Straight)

40-706-14 14,0 cm, 5 1/2"



Liston Bone Cutting Forceps (Angled)

40-707-19 20,0 cm, 8"