



Shangyi Mobile Operating Lamp YD01-4LED

1. Excellent cold light: adopting LED as light source, without temperature rising.
2. Excellent daylight quality and CRI meet the illumination needs of the doctors.
3. Long service life: LED life is up to 50,000 hours.
4. Energy saving and environment-friendly.
5. Germany OSRAM Brand LED lamp bulb with high quality lens, Each LED can be replaced individually, which guarantees maximum economy.
6. Renew Operation irradiation parts rapidly and improves doctors' working environment.
7. Illuminance equality, high shadowless rate, high definition to improve doctor's working efficiency.

YD01-4 LED Hole Type Lamp AC/DC vertical 4-reflector luminescence shadowless lamp meets the requirements for surgical operations.



Type	Shangyi Mobile Operating Lamp YD01-4LED
Illuminance(Lux)	80,000Lux
Colour temperature(k)	4800+-300
Rated power of bulb(V/W)	3.2V/1W
Service life of illuminant	50000h
Power supply voltage	AC110V-240V(50HZ-60HZ) Customize for free
Total power consumption(W)	32
Total LED bulb quantity(pcs)	32
Size of light field	160-220mm
Adjust Brightness	by control dimmer
Lamp head diameter	160mm
Illuminance depth (mm)	500mm
Emergency time	3h

Gross Weight: 54KG
Net Weight: 42KG



AKL: 21602914240