

MEDICAL GAS OUTLETS AND ADAPTERS


SPECIFICATIONS:

- Connection mode: Thread, Barb.
- Gas type: Oxygen, Air, Vacuum, Nitrous Oxide.
- 100% tested for gas leakage.
- 10,000 times Swap Fatigue Test.
- The maximum working pressure: 0,75 MPa.


FEATURES:

- Modular design.
- The main material is pure brass.
- Easy installation and well connecting.

JISS Standard



Ohmeda Standard



CIG Standard



DINN Standard



BS Standard




AGSS-Standard



JISS Adapter (Barb)



Ohmeda Adapter (Barb)



CIG Adapter (Barb)

